

Using the DEP LIMS Browser and Cost Projector

1. SCOPE

The Laboratory Information Management System (LIMS) Browser is a tool used for searching and viewing data. The Browser also can be used to create electronic or hard copy reports.

2. BROWSER ACCESS

To access the Browser, go to the *General* menu and select *LIMS Browser*.

2.1 LIMS Browser: Screen 1 – Searching for Data

Search

Search By

☐ Request ID ☐ Field ID

☒ Event ID ☐ Site Location

☐ Job ID ☐ Customer/Project

☐ Job Description ☐ Project ID

☐ Sample ID ☐ WIN/STORET #

Additional Selections

☒ None

☐ Analyses

☐ Components

☐ Analytical Group

Job Status

☒ Reported ☒ A ☐ C ☐ V ☐ X

Options

☐ Show All

☐ Show only criminal investigations

Password

Filter By Date:

☐ None

☐ Sampled

☐ Received

☐ Completed

☒ Reported

From: 12/28/2021

To: 1/27/2022

Search For...

In Database:

☒ Active ☐ Committed ☐ Legacy

Event ID

Build List

OK Cancel

Search By

Records can be searched by selecting any one of the options listed.

- Request ID: ID number generated when the event gets scheduled.
- Event ID: The event name generated when the samples are logged into the LIMS. Begins with the Customer name, e.g., WQAP.
- Job ID: Within each event there are one or more jobs. Each Job ID contains data related to a specific analytical group.
- Job Description: Refers to the description that is created in the Scheduler under "Sampling Event Description."
- Sample ID: The unique number that is assigned to a container by the LIMS system.
- Field ID: The identifier assigned to one or more containers linked to a site identifier created at the time of sample collection.
- Project ID: The LIMS Project associated with the customer for this event.
- WIN/STORET #: An identifier for the station in WIN/STORET where samples were collected. This is an optional field for LIMS and may or may not be populated at login.

If the entire search criterion number or name is not known, then a portion of it can be entered and the program will perform a wildcard search. After typing the information in the *Search By* box, click *Build List*. Select the items of interest and click OK.

Note: The search by customer/project will also bring up a project list. Click *Build List* for the projects then select one or more projects.

Additional Selections

The search can be narrowed down to specific analytes, components or analytical groups. If one of these options is selected an analyses/component/analytical group list will display after the data have been retrieved. The relevant analyses/components can then be selected.

Job Status

The data search can filter out certain job statuses. For example, if only V status jobs are needed then check only that status. For general data searches it is recommended that all statuses be checked.

Statuses:

- Reported: All jobs in the event are A status but the event may or may not have been reported to the client. This is a remnant of when the lab reported final results at the job level.

The screenshot shows the 'Search' dialog box. The 'Search By' section has 'Job ID' selected. The 'Additional Selections' section has 'Analyses' selected. The 'Job Status' section has 'Reported', 'A', 'C', 'V', and 'X' all checked. The 'Options' section has 'Show All' checked. The 'Filter By Date' section has 'Received' selected with a date range from 4/14/2003 to 4/20/2003. The 'Search For...' field contains 'TLH-2003-'. The 'In Database' section has 'Active' checked. The 'Job ID' list shows several job IDs, with 'TLH-2003-04-15-03' selected. The 'Analyses' list shows several analytes, with 'BSY-AS-MFS' selected. The 'OK' button is highlighted.


- A: The job is in status A but the event may or may not have been reported if there are other jobs in the event that have not been authorized or the report has not been certified.
- C: The analyst has uploaded the data but they have not been reviewed and authorized by a supervisor.
- V: The job does not have all the data uploaded yet.
- X: The job has been cancelled. Samples in a cancelled job do not get analyzed.

Options

- Show All: If the *Show All* box is checked, the program will retrieve all events associated with the Request IDs identified by the search criteria. This option is available when searching by Event, Job, Job Description, Sample ID, Field ID or Site Location and no additional selections (Analyses or Components) are indicated.
- Show Only Criminal Investigations: The only staff that have access to data from criminal cases are the individuals who scheduled the analyses and laboratory staff. Any others needing access may contact the Laboratory or the project manager to obtain the password. If this box is checked the search will bring up only the criminal case events.
- Filter by Date: By selecting one of these items and using a small date range you can narrow the search. For example, if you know the date collected but have no other information you can bring up a list of all samples collected on that day and then select from the list. If “None” is selected the search may encompass many years so it is best to select a type of date range when searching for data.
- In Database: There are three options: Active, Committed and Legacy. The Active database contains data for events not reported and those certified within the past 90 days. Older events are transferred to an identical set of tables called the committed database. This construct allows faster retrievals for more recent data. The Legacy database represents records from an older version of LIMS (pre-1999) and some of the search parameters are not applicable.

2.2. LIMS Browser Screen 2 – Viewing Data

After clicking OK on one or more selected items the screen below will display.




The screenshot shows the 'LIMS Browser' window. The title bar says 'LIMS Browser'. Below the title bar, there's a 'Browser' section with icons for file operations. The main area displays 'Event: CHEM-LAB-2005-12-19-01 with Status: V'. Below this, it says 'Reports Certified: None' and 'Internal PT'. Further down, it shows 'Request ID: RQ-2005-12-12-49', 'Customer: CHEM-LAB', 'Project: CHEM-QC', and 'PMAS:'. A tree view on the left shows 'Event: CHEM-LAB-2005-12-19-01' selected. The main pane displays a table of jobs:

Job ID	Status	Description	Group	Analysis Group
TLH-2005-12-19-07	V	Internal PT	CHEMISTRY	GC-Soil
TLH-2005-12-19-08	A	Internal PT	CHEMISTRY	Mercury-Water
TLH-2005-12-19-09	V	Internal PT	CHEMISTRY	Metals-Water
TLH-2005-12-19-10	V	Internal PT	CHEMISTRY	Nutrients-Liquid

At the bottom, it says 'Event: CHEM-LAB-2005-12-19-01'.

The information can be viewed by clicking on the event ID. As each new level appears the item can be clicked on to go to the next level down.

For viewing or extracting data, click on the *Report/Files* button. The screen below will appear.



The screenshot shows the 'Reports/File Formats' window. The title bar says 'Reports/File Formats'. The main area is divided into several sections:


- Report**: Contains a 'Select Event(s)' list with 'CHEM-LAB-2006-01-27-01' and 'CHEM-LAB-2006-02-17-01' selected. There is also a 'Select Job(s)' list.
- Report Title**: A text field for entering the report title.
- Detects Only**: A checkbox.
- Report Orientation**: Radio buttons for 'Portrait' and 'Landscape' (selected).
- Template**: A section for selecting fields and column order. It includes an 'Available Fields' list on the left with items like 'ANALYSIS', 'ANALYSIS DESCRIPT', 'ANALYSIS GROUP', etc. In the center, there's a 'Selected Fields/Column Order' table with columns 'Selected Field', 'Delimiter', and 'Size'. On the right, there's a 'Sort By' list.

Highlight the event(s) of interest. Then highlight the job(s) to retrieve.


Under the *Available Fields* choose the fields to include in the report. Highlight each field and click on the top arrow or double-click to move the parameter to the middle window.

The program will display the fields in the order they are selected. It will also sort by the parameters in the order they are selected unless a different sort is chosen by highlighting a parameter in the *Selected Fields* window and clicking on the top arrow by the *Sort By* window.

The size of the fields (column width) can be customized by selecting the fields of interest and clicking on “Size” in the header. Increase or decrease the number to enlarge or shrink the column width. Note that this does not affect an exported file, just the formatted PDF report.



After the fields have been selected a report can be displayed by clicking on the  *View/Print* button or choosing that menu option. The parameters can also be assigned a delimiter to create a file which can be loaded into Excel. Note that tab is the default delimiter.



To create a delimited file, highlight the parameters in the *Selected Fields* window and click on the *Delimiter* header. A window will display the options. Click the delimiter of interest.

After choosing the delimiter, click OK. Click on the  *Create Delimited File* button or choose the menu

option. A dialog box will be displayed, prompting the user to save the file. Select the directory for saving, name the file and click OK to save.

The customized field selections can be saved as a template. After the fields, sort  by and delimiter selections have been made, click on the icon pictured to the left. In the window that pops up name the template.

To use an existing template, highlight the events/jobs of interest and click  the *Open Template* button. Select the template of interest and click OK.

Other features

To view results of only values above the detection limits then check the box *Detects Only*.

To exclude results that have been cancelled, check the *Exclude Cancelled* box.

To retrieve field parameters that were entered using the WIN Data Merge Tool, check the *Additional Field Parameters* box.

Viewing Final Reports

For events that have been reported you can view a copy of the signed final report using the LIMS browser.

To view the final report, highlight the event of interest. Click on *Report* in the upper left corner. On the drop-down menu that appears click on *Show Final Report*.

If there are both chemistry and biology analyses in the event there will be two separate reports. The selection box to choose chemistry or biology will appear if both reports are available. A PDF file of the final report will open.

3.0 LIMS COST PROJECTOR

The Cost Projector can be used to obtain the value of scheduled samples (samples requested but not yet submitted) and events (samples that have been submitted to the lab).

3.1 From *General* menu, select *Cost Projector*.

3.2 To get the cost of samples which have been submitted to the lab:

Click on the Event button. Select or enter the event name. If you don't know the event name but know the request ID (RQ), use the Browser to look up the event name.

When using the RQ, you may get an incorrect cost by using the Scheduler button. Cost by RQ gives the cost of what was scheduled not what was submitted to the lab.

- 3.3 To get the cost of samples that have been scheduled but not submitted to the lab:

Click on the Scheduler button and enter or select the RQ of interest.

Appendix of Changes

- 02/11/2016 Updated language and substituted screen shots of buttons to replace verbal descriptions.
- 02/14/2017 Added language and an updated screen shot to describe the search by WIN/STORET #.
- 01/27/2022 Updated screen shot to display the search by Analytical Group option.