Effective Date: 04/22/2026
Revision Date: 04/22/2026
Revision Author: A. Terrane

Balance Daily Use Record Mettler XPR105DR, Micro Lab
(Weights Bracketing Measured Weight)

TA-MICRO-XPR-LOG-1.0

Check Weight IDs
BIO-R-WEIGHTS-02 = R02
BIO-R-WEIGHTS-04 = R04
BIO-R-WEIGHTS-01 = R05
BIO-S-WEIGHTS-01 = S01

STATE of FLORIDA
DEPARTMENT of ENVIRONMENTAL PROTECTION
BIOLOGY SECTION

Mettler XPR105DR serial # C438788153

Annual Calibration Required

Date/Time	Material Weight or Weight Range	Lower Weight Value	Lower Weight Measured Value	Upper Weight Value	Upper Weight Measured Value	Check Weight ID	Weighed by: (Initials)	Comments

Load limits: 1-40g (high precision, 0.00001g readability)
41-120g (standard precision, 0.0001g readability)
DO NOT EXCEED 120g TOTAL LOAD

NOTE: Record g for grams; mg for milligrams
One lower value and one upper value check-weight can be used to bracket an event consisting of numerous measurements.
Additional check-weight measurements are required (per day) if more than two hours elapse between measurements.

Accuracy: Approx. ±0.00002 g at ~1 g; up to ±0.00083 g at 120 g (accuracy decreases as load increases)

Comments: