Preparation of Lugol's Solution

(Modified from Standard Methods 10200B.2)

1. SCOPE AND APPLICATION

1.1 This method details the preparation of Lugol's solution, which is used to preserve phytoplankton samples.

1.2 This preservative is used for samples with LIMS test codes PKD-QN-C and DTY-QN-C.

2. SUMMARY OF METHOD

2.1 Deionized (DI) water, acetic acid, potassium iodide, and iodine crystals are mixed in the proper proportions, labeled, and stored.

3. APPARATUS AND EQUIPMENT

3.1 Mettler Toledo balance, PB403-S/FACT or PC440DeltaRange®
3.2 Weigh boat
3.3 Spatula
3.4 Magnetic stir plate
3.5 Erlenmeyer flask, 500-mL
3.6 Teflon stir bar
3.7 Glass amber bottle, 500-mL

4. REAGENTS AND CHEMICALS: AMOUNTS SHOWN WILL MAKE 450 ML OF LUGOL'S SOLUTION.

4.1 Potassium iodide (KI), 40 g
4.2 Iodine crystals, 20 g
4.3 Glacial acetic acid, 40 mL
4.4 Deionized (DI) water, 400 mL

5. METHODS

5.1 Work in a fume hood. In the 500-mL Erlenmeyer flask, add 40 mL of acetic acid to 400 mL of DI water.

5.2 Place the 500-mL Erlenmeyer flask on a magnetic stir plate and put the stir bar inside of flask.

5.3 Add 40 g of potassium iodide and 20 g of iodine crystals to the acid/water solution.

5.4 Stir solution for 30 minutes. Check to make sure all of the crystals have dissolved.

Note: If needed, the solution can be stirred for up to an hour.

5.5 Pour mixed solution into 500-mL glass amber bottle. Label container with "Lugol's Soln", the prep date, analyst initials, and a batch number as generated per step 5.6.

5.6 Generate a batch number in the Standard Preparation Tracker module in LIMS. Record the batch number (by hand or print a sticker) on the container.

The OFFICIAL COPY is the on-line version. All other copies are considered unofficial and uncontrolled.

5.6.1 To create a batch number: LIMS --> Laboratory --> Standard Preparation Tracker: Select Group = Bio-Taxonomy. Select Standard = New or Duplicate. If using New, enter all the information required. If using Duplicate, select the Identity of the solution type you are duplicating. Review and complete the information required for accuracy. A serial number (batch number) will be generated. Record this number. From Tools, print a label (do not print from Standard if prompted because the sticker will not include all the necessary information).

Note: Be sure and record the proper amounts of reagents to avoid confusion.

5.7 Store the container under a fume hood in the lab.

6. QUALITY CONTROL AND DOCUMENTATION

6.1 LIMS Standard Preparation Tracker module

6.2 To correct raw data entries, place a single line through the incorrect entry, write the corrected entry near the error and initial the correction. If the correction requires an explanation, write a comment and reference the error.

6.3 See the FDEP Laboratory Quality Manual

7. SAFETY

7.1 Lab coat, latex or vinyl gloves, and eye protection are required when working with chemicals and unprocessed samples.

7.2 Avoid breathing fumes or exposing co-workers to them. Use proper safety equipment, a lab hood or chemical ventilation system, and keep face away from open containers. Work in a well ventilated area.

7.3 Reagent bottles must be labeled with the name and concentration of the reagent.

7.4 Refer to the Material Safety Data Sheet (MSDS) for information about reagents and chemicals in this SOP.

7.5 Refer to the FDEP Laboratory Health and Safety Plan.

7.6 Refer to the FDEP Laboratory Quality Manual.

8. POLLUTION PREVENTION AND WASTE MANAGEMENT

8.1 Refer to the FDEP Laboratory Health and Safety Plan.

9. REFERNECES

9.1 FDEP Laboratory Quality Manual, http://www.dep.state.fl.us/labs/library/lab_sops.htm

9.2 FDEP Laboratory Health and Safety Plan, https://www.fldepnet.org/content/dear/bureau-laboratories-safety-documents

Appendix of Changes

07/12/2010 Reformatted the layout following SOP LB-001.

06/27/2013 Revised step 5.6 to use Standard Preparation Tracker module in LIMS.

The OFFICIAL COPY is the on-line version. All other copies are considered unofficial and uncontrolled.

07/01/2015 Formatted for 508 compliance. Added References section.

08/01/2016 Removed reference to lugol's logbook; using LIMS Prep Tracker to generate batch number.

The OFFICIAL COPY is the on-line version. All other copies are considered unofficial and uncontrolled.