

Certificate of Analysis List

For request number 1085788

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description	
2660501 NTU	7207	2659405	7223	STABLCAL [®] FORMAZIN STANDARD	800
2684701	7163	2659405	7223	<0.1 NTU Calibration Solution	
2684801 NTU	7188	2659405	7223	STABLCAL [®] FORMAZIN STANDARD	20
2684901 NTU	7202	2659405	7223	STABLCAL [®] FORMAZIN STANDARD	100

Total Enclosures: 4

***Certificate of Analysis******This is a Component of 2659405 lot A7223***

Page 1

COMMODITY: **STABLCAL® FORMAZIN STANDARD 800 NTU**COMMODITY NUMBER: **2660501**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A7207****8/3/2017****8/8/2017**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	760 to 840 NTU	826.7 NTU

The expiration date is Nov 2018

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

***Certificate of Analysis******This is a Component of 2659405 lot A7223***

Page 1

COMMODITY: <0.1 NTU Calibration Solution

COMMODITY NUMBER: 2684701

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: A7163

6/16/2017

6/20/2017

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	0 to 0.1 NTU	0.06 NTU

The expiration date is Sep 2018

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

***Certificate of Analysis******This is a Component of 2659405 lot A7223***

Page 1

COMMODITY: **STABLCAL® FORMAZIN STANDARD 20 NTU**COMMODITY NUMBER: **2684801**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A7188****7/10/2017****7/27/2017**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	19 to 21 NTU	20.4 NTU

The expiration date is Oct 2018

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

***Certificate of Analysis******This is a Component of 2659405 lot A7223***

Page 1

COMMODITY: **STABLCAL® FORMAZIN STANDARD 100 NTU**COMMODITY NUMBER: **2684901**

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: **A7202****7/31/2017****8/4/2017**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	95 to 105 NTU	102.4 NTU

The expiration date is Oct 2018

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist