

Certificate of Analysis List

For request number 1889438

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2684701	3144	N/A	N/A	<0.1 NTU Calibration Solution

Total Enclosures: 1

***Certificate of Analysis***

Page 1

COMMODITY: <0.1 NTU Calibration Solution

COMMODITY NUMBER: 2684701

MANUFACTURE DATE:

DATE OF ANALYSIS:

LOT NUMBER: A3144

6/1/2023

6/5/2023

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	0 to 0.1 NTU	0.04 NTU

The expiration date is Sep 2024

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

Certificate of Analysis List

For request number 1889441

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2684801 STANDARD	3158 20 NTU	N/A	N/A	STABLCAL sup TS sup0 FORMAZIN

Total Enclosures: 1

***Certificate of Analysis***

Page 1

COMMODITY: **STABLCAL|sup|TS|sup0 FORMAZIN STANDARD 20 NTU**
COMMODITY NUMBER: **2684801** MANUFACTURE DATE: DATE OF ANALYSIS:
LOT NUMBER: **A3158** **6/12/2023** **6/20/2023**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	19 to 21 NTU	20.4 NTU

The expiration date is Sep 2024

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

Certificate of Analysis List

For request number 1889442

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2684901 STANDARD	3157 100 NTU	N/A	N/A	STABLCAL sup TS sup0 FORMAZIN

Total Enclosures: 1

*Certificate of Analysis*

Page 1

COMMODITY: **STABLCAL|sup|TS|sup0 FORMAZIN STANDARD** 100 NTU
COMMODITY NUMBER: **2684901** MANUFACTURE DATE: DATE OF ANALYSIS:
LOT NUMBER: **A3157** 6/7/2023 6/16/2023

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	95 to 105 NTU	101.0 NTU

The expiration date is Sep 2024

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist

Certificate of Analysis List

For request number 1889434

Catalog Number Entered	Lot Number Entered	Related Catalog Number	Related Lot Code	Description
2660501 STANDARD 800 NTU	3157	N/A	N/A	STABLCAL sup TS sup0 FORMAZIN

Total Enclosures: 1

***Certificate of Analysis***

Page 1

COMMODITY: **STABLCAL|sup|TS|sup0 FORMAZIN STANDARD 800 NTU**
COMMODITY NUMBER: **2660501** MANUFACTURE DATE: DATE OF ANALYSIS:
LOT NUMBER: **A3157** **6/8/2023** **6/20/2023**

<i>TEST</i>	<i>SPECIFICATIONS</i>	<i>RESULTS</i>
Turbidity	760 to 840 NTU	824.0 NTU

The expiration date is Sep 2024

Formazin and StablCal® solutions provided by Hach are not NIST traceable because the NIST does not carry turbidity standards. However, the use of Formazin and StablCal® as used in Hach method 8195 are accepted by the EPA as a primary standard to be used in the calibration of turbidity instruments.

Certified by _____

Scott Als
Analytical Services Chemist