

Chemical Analysis Report

NE-DIST-2017-12-20-03

Florida Department of Environmental Protection
Central Laboratory
2600 Blair Stone Road
Tallahassee, FL 32399-2400
DOH Accreditation E31780

Event Description: **Suwanee spring SCI**
Request ID: **RQ-2017-12-18-06**
Customer: **NE-DIST**
Project ID: **SMP**

Send Reports to:
FL Dept. of Environmental Protection
8800 Baymeadows Way West
Jacksonville, FL 32256-7577
Attn: Jeremy Parrish

For additional information please contact
Colin Wright, Ph.D.
Liang-Tsair Lin, Ph.D.
Thekkekalathil Chandrasekhar, Ph.D, QA Officer
Phone (850) 245-8085

Certified by: Thekkekalathil Chandrasekhar, Chemist Administrator

Date Certified: 05-JAN-2018 11:15



Case Narrative

Unless otherwise noted, all samples included in this report were received in accordance with protocols referenced in Chapter 62-160, Florida Administrative Code (F.A.C.). Results published in this report pertain only to the samples as submitted to, and received by the laboratory. All times in this report are adjusted to the applicable Eastern Time Zone (EST or EDT).

Results for the following analytical groups are included in this report: Nutrients and Pesticides.

Scientific notation may be used in reporting very large or small values. Values reported using scientific notation will take the form of the following example: 1.3E+03, which is equivalent to 1.3×10^3 or 1300.

Unless otherwise noted, analytical values for soil and sediment samples are reported on a dry weight basis, and analytical values for waste and tissue samples are reported on a wet weight basis.

Results for TNI accredited tests met requirements established by The NELAC Institute. A double asterisk (**) is used to indicate an analyte/matrix/method for which the laboratory is not TNI accredited or where accreditation for that field of testing is not applicable.

Any significant anomalies or deviations from established protocols are documented in Non-Conformance Reports, which, where appropriate, are included within this analytical report. Additional comments related to specific analytical tests may be included as remarks following the analytical results for each sample. Such comments and remarks are for informational purposes only and are not intended to convey judgement about the usability of the reported data.

A quality control report on the performance of the test method for the submitted samples is included. Uncertainty associated with the analytical results contained in this report can be estimated from the reported quality assurance results and from published quality control acceptance limits for each analytical test. Matrix quality control results (matrix spike recoveries and matrix sample precision) pertain only to the matrix sample tested and do not necessarily reflect test method performance for other samples.

Typical matrix quality control (QC) measurements may include matrix spike recovery, matrix spike duplicate recovery, matrix spike precision and matrix sample precision. Not all matrix QC results may be available or reportable; where they are not an explanation is provided. Typical reasons for unavailable QC results include, but are not limited to, a) insufficient matrix sample to perform some or all QC measurements; b) analyte concentration in the sample replicated was too low for a meaningful measurement of precision and c) analyte concentration in the matrix sample spiked was too high (relative to the amount of analyte spiked) for a meaningful measurement of recovery. Where matrix QC results are unavailable, other method performance metrics (e.g., LCS recovery, LCS precision, surrogate recovery) may be used to assess performance of the method. Comments explaining any missing QC measurements are not intended to convey any adverse conclusions about the quality of the reported data.

Precision is reported as relative percent difference unless otherwise noted.

Quality Control codes as defined below may be used in this report to indicate results that are associated with one or more quality control elements which did not fall within established test method criteria. Such results may be qualified as estimates using a J qualifier as required by 62-160 F.A.C. Explanations are included in the report for any results that were reported as estimates for other reasons.

QC Codes used in this report may include:

- LCS – Recovery for the batch Laboratory Control Sample (LCS) was outside existing control limits;
- MS – Recovery for the batch matrix spike (MS) was outside existing control limits;
- CCV – Recovery for a continuing calibration verification (CCV) standard was outside existing control limits;
- SUR – Recovery of a surrogate (SUR) for associated analytes was outside existing control limits;
- RPD – The precision, measured as relative percent difference (RPD), of batch replicate measurements was outside existing control limits;
- RSD – The precision, measured as relative standard deviation (RSD), of batch replicate measurements was outside existing control limits;
- SMP – Sample - used precision derived from replicate analyses of a sample;

The following data qualifiers are used, where applicable, in this report as specified in 62-160 F.A.C.

- A - Value reported is the mean of two or more determinations.
- B - Results based on colony counts outside the acceptable range.
- I - The reported value is between the laboratory method detection limit and the laboratory practical quantitation limit.
- J - Estimated value and/or the analysis did not meet established quality control criteria.
- K - Actual value is known to be less than value given.
- L - Actual value is known to be greater than value given.
- N - Presumptive evidence of presence of material.
- O - Sampled, but analysis lost or not performed.
- Q - Sample held beyond normal holding time.
- T - Value reported is less than the criterion of detection.
- U - Material was analyzed for but not detected. The reported value is the method detection limit for the sample analyzed.
- V - Analyte was detected in both sample and method blank.
- X - Too few individuals to calculate SCI value.
- Y - The laboratory analysis was from an unpreserved or improperly preserved sample. The data may not be accurate.
- Z - Colonies were too numerous to count (TNTC).

Quality control information from overflow laboratories may not be included in this report. Please refer to the associated report from the overflow laboratory for additional information.

Sample Location: WACASSASSA R @ SR 24

Collection Date/Time: 12/19/2017 11:20

Field ID: G1NE0829 12192017

Matrix: W-SURF-FRH

Sample ID	Ref. Method	Component	Result	Code	Units	Batch ID	QC Codes
1955404	EPA 180.1 Rev. 2.0	Turbidity	0.75		NTU	P337305	
	EPA 300.0 Rev. 2.1	Chloride	8.5		mg Cl/L	P337392	
		Sulfate	1.4		mg SO4/L	P337396	
	SM 2120 B	Color (true)	180		PCU	P337292	
	SM 2320 B-97	Total Alkalinity	79		mg CaCO3/L	P337489	
	SM 2540 C-97	TDS	138	A	mg/L	P337670	
	SM 2540 D-97	TSS	2	IA	mg/L	P337674	
	SM 4500 F-C-97	Fluoride	0.059	I	mg F/L	P337492	
1955405	EPA 350.1 Rev. 2.0	Ammonia-N	0.006		mg N/L	P337700	
	EPA 351.2 Rev. 2.0	Kjeldahl Nitrogen	0.68		mg N/L	P337403	
	EPA 353.2 Rev. 2.0	NO2NO3-N	0.004	U	mg N/L	P337407	
	EPA 365.1 Rev. 2.0	Total-P	0.026		mg P/L	P337417	
	SM 5310 B-00	Organic Carbon	19		mg C/L	P337473	
1955406	EPA 365.1 Rev. 2.0 dissolved	O-Phosphate-P	0.014		mg P/L	P337299	

Ref. Method and Comment:

SM 2540 D-97: Batch relative percent difference is unavailable because of low analyte concentration in the QC sample.

Sample Location: WACASSASSA R @ SR 24

Collection Date/Time: 12/19/2017 11:20

Field ID: G1NE0829 12192017

Matrix: W-SURF-FRH

Sample ID	Ref. Method	Component	Result	Code	Units	Batch ID	QC Codes
1955410	EPA 8321B	Sucralose**	0.010	U	ug/L	P337235	
	USGS O-2060-01	2,4-D	0.0020	U	ug/L	P337256	
		Diuron	8.0E-04	U	ug/L	P337256	
		Fenuron	0.0080	U	ug/L	P337256	
		Imidacloprid	8.0E-04	U	ug/L	P337256	MS
		Bentazon	8.0E-04	U	ug/L	P337256	
		Linuron	0.0040	U	ug/L	P337256	
		Acetaminophen**	0.0080	U	ug/L	P337256	
		MCPD	0.0020	U	ug/L	P337256	
		Carbamazepine**	4.0E-04	U	ug/L	P337256	
		Triclopyr**	0.0040	U	ug/L	P337256	
		Primidone**	0.0040	U	ug/L	P337256	
		Fluridone**	4.0E-04	U	ug/L	P337256	

Ref. Method and Comment:

USGS O-2060-01: MS accuracy and precision for bentazon could not be assessed due to high concentration of parameter in the spiked sample.

Quality Assurance Report

Method Blank Results

Reference Method: EPA 180.1 Rev. 2.0

Batch ID: P337305

Component	Result	Code	Units
Turbidity	0.10	U	NTU

Quality Assurance Report Method Blank Results

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P337392

Component	Result	Code	Units
Chloride	0.20	U	mg Cl/L

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P337396

Component	Result	Code	Units
Sulfate	0.20	U	mg SO4/L

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P337700

Component	Result	Code	Units
Ammonia-N	0.002	U	mg N/L

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P337403

Component	Result	Code	Units
Kjeldahl Nitrogen	0.080	U	mg N/L

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P337407

Component	Result	Code	Units
NO2NO3-N	0.004	U	mg N/L

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P337417

Component	Result	Code	Units
Total-P	0.002	U	mg P/L

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P337299

Component	Result	Code	Units
O-Phosphate-P	0.004	U	mg P/L

Reference Method: EPA 8321B
Batch ID: P337235

Component	Result	Code	Units
Sucralose	0.010	U	ug/L

Reference Method: SM 2120 B
Batch ID: P337292

Component	Result	Code	Units
Color (true)	2.5	U	PCU

Quality Assurance Report Method Blank Results

Reference Method: SM 2320 B-97
 Batch ID: P337489

Component	Result	Code	Units
Total Alkalinity	0.65	U	mg CaCO ₃ /L

Reference Method: SM 2540 C-97
 Batch ID: P337670

Component	Result	Code	Units
TDS	15	U	mg/L

Reference Method: SM 2540 D-97
 Batch ID: P337674

Component	Result	Code	Units
TSS	2	U	mg/L

Reference Method: SM 4500 F-C-97
 Batch ID: P337492

Component	Result	Code	Units
Fluoride	0.035	U	mg F/L

Reference Method: SM 5310 B-00
 Batch ID: P337473

Component	Result	Code	Units
Organic Carbon	0.50	U	mg C/L

Reference Method: USGS O-2060-01
 Batch ID: P337256

Component	Result	Code	Units
2,4-D	0.0020	U	ug/L
Acetaminophen	0.0080	U	ug/L
Bentazon	8.0E-04	U	ug/L
Carbamazepine	4.0E-04	U	ug/L
Diuron	8.0E-04	U	ug/L
Fenuron	0.0080	U	ug/L
Fluridone	4.0E-04	U	ug/L
Imidacloprid	8.0E-04	U	ug/L
Linuron	0.0040	U	ug/L
MCP	0.0020	U	ug/L
Primidone	0.0040	U	ug/L
Triclopyr	0.0040	U	ug/L

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: EPA 300.0 Rev. 2.1
 Batch ID: P337392

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Chloride	105		P	90 - 110

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P337396

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Sulfate	100		P	90 - 110

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P337700

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Ammonia-N	103		P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P337403

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Kjeldahl Nitrogen	99.9		P	90 - 110

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P337407

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
NO2NO3-N	102		P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P337417

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Total-P	101		P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P337299

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
O-Phosphate-P	99.4		P	90 - 110

Reference Method: EPA 8321B
Batch ID: P337235

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Sucralose	89.5		P	40 - 150

Reference Method: SM 2320 B-97
Batch ID: P337489

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Total Alkalinity	102		P	95 - 108

Reference Method: SM 4500 F-C-97
Batch ID: P337492

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Fluoride	94.2		P	92 - 103

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: SM 5310 B-00
 Batch ID: P337473

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Organic Carbon	95.2		P	85 - 115

Reference Method: USGS O-2060-01
 Batch ID: P337256

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
2,4-D	97.3		P	40 - 150
Acetaminophen	100		P	30 - 150
Bentazon	89.3		P	30 - 150
Carbamazepine	101		P	40 - 150
Diuron	90.8		P	40 - 150
Fenuron	113		P	40 - 150
Fluridone	109		P	40 - 150
Imidacloprid	105		P	40 - 160
Linuron	99.6		P	40 - 150
MCP	89.3		P	40 - 150
Primidone	103		P	40 - 150
Triclopyr	88.2		P	40 - 150

Quality Assurance Report Matrix Spike Accuracy

Reference Method: EPA 300.0 Rev. 2.1
 Batch ID: P337392

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955267	Chloride	105		P	90 - 110
1955302	Chloride	103		P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1
 Batch ID: P337396

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955267	Sulfate	101		P	90 - 110
1955302	Sulfate	101		P	90 - 110

Reference Method: EPA 350.1 Rev. 2.0
 Batch ID: P337700

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955339	Ammonia-N	101	102	P/P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0
 Batch ID: P337403

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955339	Kjeldahl Nitrogen	104	98.4	P/P	90 - 110

Quality Assurance Report Matrix Spike Accuracy

Reference Method: EPA 353.2 Rev. 2.0
 Batch ID: P337407

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955317	NO2NO3-N	104	104	P/P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0
 Batch ID: P337417

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955532	Total-P	96.2	100	P/P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0 dissolved
 Batch ID: P337299

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955306	O-Phosphate-P	99.0	98.9	P/P	90 - 110

Reference Method: EPA 8321B
 Batch ID: P337235

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955251	Sucralose	105	108	P/P	40 - 150

Reference Method: SM 4500 F-C-97
 Batch ID: P337492

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955487	Fluoride	94.9	96.1	P/P	91 - 105

Reference Method: SM 5310 B-00
 Batch ID: P337473

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955140	Organic Carbon	97.1	99.4	P/P	85 - 115

Reference Method: USGS O-2060-01
 Batch ID: P337256

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1955250	2,4-D	113	116	P/P	40 - 150
1955250	Acetaminophen	115	109	P/P	30 - 150
1955250	Carbamazepine	112	104	P/P	40 - 150
1955250	Diuron	94.9	89.2	P/P	40 - 150
1955250	Fenuron	124	93.5	P/P	40 - 150
1955250	Fluridone	96.2	88.2	P/P	40 - 150
1955250	Imidacloprid	172	150	F/P	40 - 160
1955250	Linuron	101	89.3	P/P	40 - 150
1955250	MCP	86.6	86.8	P/P	40 - 150
1955250	Primidone	114	106	P/P	40 - 150
1955250	Triclopyr	91.5	88.4	P/P	40 - 150

Quality Assurance Report Precision

Quality Assurance Report Precision

Reference Method: EPA 180.1 Rev. 2.0
Batch ID: P337305

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955349	Turbidity	0.0	Sample	P	0 - 20

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P337392

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955302	Chloride	1.99	Sample	P	0 - 20

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P337396

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955302	Sulfate	2.91	Sample	P	0 - 20

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P337700

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955339	Ammonia-N	1.22	Spike	P	0 - 20

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P337403

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955339	Kjeldahl Nitrogen	4.00	Spike	P	0 - 20

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P337407

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955317	NO2NO3-N	0.0	Spike	P	0 - 20

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P337417

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955532	Total-P	3.81	Spike	P	0 - 20

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P337299

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955306	O-Phosphate-P	0.0969	Spike	P	0 - 20

Quality Assurance Report Precision

Reference Method: EPA 8321B
Batch ID: P337235

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955251	Sucralose	1.64	Spike	P	0 - 30

Reference Method: SM 2120 B
Batch ID: P337292

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955341	Color (true)	3.93	Sample	P	0 - 20

Reference Method: SM 2320 B-97
Batch ID: P337489

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955487	Total Alkalinity	0.353	Sample	P	0 - 4.8

Reference Method: SM 2540 C-97
Batch ID: P337670

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955404	TDS	0.727	Sample	P	0 - 10

Reference Method: SM 4500 F-C-97
Batch ID: P337492

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955487	Fluoride	0.953	Spike	P	0 - 3.0

Reference Method: SM 5310 B-00
Batch ID: P337473

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955140	Organic Carbon	2.27	Spike	P	0 - 20

Reference Method: USGS O-2060-01
Batch ID: P337256

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955250	2,4-D	1.97	Spike	P	0 - 30
1955250	Acetaminophen	5.36	Spike	P	0 - 30
1955250	Carbamazepine	6.52	Spike	P	0 - 30
1955250	Diuron	5.88	Spike	P	0 - 30
1955250	Fenuron	27.9	Spike	P	0 - 30
1955250	Fluridone	6.18	Spike	P	0 - 30
1955250	Imidacloprid	12.5	Spike	P	0 - 30
1955250	Linuron	12.4	Spike	P	0 - 30
1955250	MCP	0.138	Spike	P	0 - 30

Quality Assurance Report Precision

Reference Method: USGS O-2060-01
Batch ID: P337256

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1955250	Primidone	7.44	Spike	P	0 - 30
1955250	Triclopyr	3.39	Spike	P	0 - 30

* Sample, spike and/or laboratory control sample precision (LCS) is reported.
Replicate spike precision may be reported when sample results are below quantifiable levels.

Quality Assurance Report Surrogates

Lab Sample ID: 1955410
Field Sample ID: G1NE0829 12192017

Reference Method	Surrogate	% Rec.	Pass/Fail	Control Limits
EPA 8321B	Acetaminophen-d4	120	P	30 - 160
EPA 8321B	Carbamazepine-d10	108	P	30 - 160
EPA 8321B	Diuron-d6	92.3	P	30 - 160
EPA 8321B	Primidone-d5	101	P	30 - 160
EPA 8321B	Sucralose-d6	126	P	40 - 160
USGS O-2060-01	2,4-D-d3	103	P	30 - 160
USGS O-2060-01	Acetaminophen-d4	120	P	30 - 160
USGS O-2060-01	Carbamazepine-d10	108	P	30 - 160
USGS O-2060-01	Diuron-d6	92.3	P	30 - 160
USGS O-2060-01	Primidone-d5	101	P	30 - 160
USGS O-2060-01	Sucralose-d6	126	P	40 - 160

Quality Assurance Report Calibration Verification

Reference Method: SM 2120 B
Run ID: A81043
Included Lab Sample IDs: 1955404

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Color (true)	99.5	99.7	P/P	85 - 115

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Run ID: A81046
Included Lab Sample IDs: 1955406

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
O-Phosphate-P	98.3	98.7	P/P	90 - 110

Reference Method: EPA 180.1 Rev. 2.0
Run ID: A81048
Included Lab Sample IDs: 1955404

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Turbidity	100	100	P/P	95 - 105

Quality Assurance Report Calibration Verification

Reference Method: EPA 300.0 Rev. 2.1

Run ID: A81087

Included Lab Sample IDs: 1955404

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Chloride	107	105	P/P	90 - 110
Sulfate	103	102	P/P	90 - 110

Reference Method: EPA 353.2 Rev. 2.0

Run ID: A81092

Included Lab Sample IDs: 1955405

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
NO2NO3-N	101	100	P/P	90 - 110

Reference Method: USGS O-2060-01

Run ID: A81093

Included Lab Sample IDs: 1955410

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
2,4-D	96.9	99.8	P/P	50 - 150
Acetaminophen	105	100	P/P	60 - 150
Bentazon	99.3	97.6	P/P	50 - 150
Carbamazepine	103	103	P/P	60 - 150
Diuron	102	99.7	P/P	60 - 150
Fenuron	104	102	P/P	60 - 150
Fluridone	105	105	P/P	60 - 160
Imidacloprid	96.0	96.8	P/P	60 - 150
Linuron	102	101	P/P	60 - 150
MCPP	95.8	91.0	P/P	50 - 150
Primidone	103	103	P/P	60 - 150
Triclopyr	97.2	85.8	P/P	50 - 150

Reference Method: EPA 365.1 Rev. 2.0

Run ID: A81121

Included Lab Sample IDs: 1955405

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Total-P	101	102	P/P	90 - 110

Reference Method: SM 5310 B-00

Run ID: A81124

Included Lab Sample IDs: 1955405

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Organic Carbon	96.7	97.0	P/P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0

Run ID: A81127

Included Lab Sample IDs: 1955405

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Kjeldahl Nitrogen	105	107	P/P	90 - 110

Quality Assurance Report Calibration Verification

Reference Method: SM 2320 B-97
 Run ID: A81134
 Included Lab Sample IDs: 1955404

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Total Alkalinity	101	98.5	P/P	90 - 110

Reference Method: SM 4500 F-C-97
 Run ID: A81134
 Included Lab Sample IDs: 1955404

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Fluoride	98.6	100	P/P	90 - 110

Reference Method: EPA 8321B
 Run ID: A81137
 Included Lab Sample IDs: 1955410

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Sucralose	88.6	74.6	P/P	60 - 150

Reference Method: EPA 350.1 Rev. 2.0
 Run ID: A81225
 Included Lab Sample IDs: 1955405

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Ammonia-N	102	102	P/P	90 - 110

* Pass/Fail determinations are made for each bracketing calibration verification check.
 Control limits for initial calibration checks may be different from those for continuing checks, depending on method requirements.
 Where they are different, both control limits are provided.

Quality Assurance Report Summary

Ref. Method	Analyte	LCS % Recovery	MS % Recovery	LCS	Precision SMP	MS
EPA 180.1 Rev. 2.0	Turbidity				0.0	
EPA 300.0 Rev. 2.1	Chloride	105	103	105	1.99	
	Sulfate	100	101	101	2.91	
EPA 350.1 Rev. 2.0	Ammonia-N	103	101	102		1.22
EPA 351.2 Rev. 2.0	Kjeldahl Nitrogen	99.9	104	98.4		4.00
EPA 353.2 Rev. 2.0	NO2NO3-N	102	104	104		0.0
EPA 365.1 Rev. 2.0	Total-P	101	96.2	100		3.81
EPA 365.1 Rev. 2.0 dissolved	O-Phosphate-P	99.4	99.0	98.9		0.0969
EPA 8321B	Sucralose	89.5	105	108		1.64
SM 2120 B	Color (true)				3.93	
SM 2320 B-97	Total Alkalinity	102			0.353	
SM 2540 C-97	TDS				0.727	
SM 4500 F-C-97	Fluoride	94.2	94.9	96.1		0.953
SM 5310 B-00	Organic Carbon	95.2	97.1	99.4		2.27
USGS O-2060-01	2,4-D	97.3	113	116		1.97
	Acetaminophen	100	115	109		5.36
	Bentazon	89.3				
	Carbamazepine	101	112	104		6.52

Quality Assurance Report Summary

Ref. Method	Analyte	LCS % Recovery	MS % Recovery		Precision SMP	MS
			LCS			
USGS O-2060-01	Diuron	90.8	94.9	89.2		5.88
	Fenuron	113	124	93.5		27.9
	Fluridone	109	96.2	88.2		6.18
	Imidacloprid	105	172	150		12.5
	Linuron	99.6	101	89.3		12.4
	MCP	89.3	86.6	86.8		0.138
	Primidone	103	114	106		7.44
	Triclopyr	88.2	91.5	88.4		3.39

Reference Method Descriptions

Method / DoH Cert #	Description	Associated Samples
EPA 180.1 Rev. 2.0 / E31780	Turbidity in aqueous matrices	1955404
EPA 300.0 Rev. 2.1 / E31780	Chloride in aqueous matrices	1955404
EPA 300.0 Rev. 2.1 / E31780	Sulfate in aqueous matrices	1955404
EPA 350.1 Rev. 2.0 / E31780	Ammonia in aqueous matrices as mg N/L	1955405
EPA 351.2 Rev. 2.0 / E31780	Total Kjeldahl Nitrogen in aqueous matrices	1955405
EPA 353.2 Rev. 2.0 / E31780	Nitrite/Nitrate in aqueous matrices as mg N/L	1955405
EPA 365.1 Rev. 2.0 / E31780	Total Phosphorus in aqueous matrices as mg P/L	1955405
EPA 365.1 Rev. 2.0 dissolved / E31780	Ortho-phosphate, dissolved, in filtered, aqueous matrices as mg P/L	1955406
EPA 8321B / E31780	Analysis of sucralose in water matrices by HPLC/MS/MS	1955410
SM 2120 B / E31780	True Color measured at 450 nm	1955404
SM 2320 B-97 / E31780	Alkalinity in aqueous matrices as mg CaCO ₃ /L	1955404
SM 2540 C-97 / E31780	Total Dissolved Solids in aqueous matrices	1955404
SM 2540 D-97 / E31780	Total Suspended Solids in aqueous matrices	1955404
SM 4500 F-C-97 / E31780	Fluoride in aqueous matrices, without distillation (not valid for NPDES monitoring)	1955404
SM 5310 B-00 / E31780	Nonpurgeable Organic Carbon in aqueous matrices	1955405
USGS O-2060-01 / E31780	Chlorinated (phenoxy acid) herbicides in water matrices by HPLC/MS/MS	1955410
USGS O-2060-01 / E31780	Urea herbicides and other chemicals in water matrices by HPLC/MS/MS	1955410

Preparation and Analysis Log

Ref. Method	Received Date	Prep Date/Time	Prepared By	Analysis Date/Time	Analyzed By	Associated Samples
EPA 180.1 Rev. 2.0	12/20/2017			12/20/2017 16:23	DeAsia Armster	1955404
EPA 300.0 Rev. 2.1	12/20/2017			12/21/2017	Julia Storbeck	1955404
EPA 350.1 Rev. 2.0	12/20/2017			12/29/2017	Sofia Elnemr	1955405
EPA 351.2 Rev. 2.0	12/20/2017	12/22/2017	Adrian Shelby	12/26/2017	Joseph Suarez	1955405
EPA 353.2 Rev. 2.0	12/20/2017			12/22/2017	Lipi Saha	1955405
EPA 365.1 Rev. 2.0	12/20/2017	12/22/2017	Sima Feizi	12/26/2017	Lipi Saha	1955405
EPA 365.1 Rev. 2.0 dissolved	12/20/2017			12/20/2017 15:54	Megan Treadwell	1955406
EPA 8321B	12/20/2017	12/20/2017	Juyoung Kim	12/26/2017	Juyoung Kim	1955410
SM 2120 B	12/20/2017			12/20/2017 15:02	Tanya Welborn	1955404
SM 2320 B-97	12/20/2017			12/27/2017	Dale L. Simmons	1955404
SM 2540 C-97	12/20/2017			12/25/2017	Yijie Li	1955404
SM 2540 D-97	12/20/2017			12/25/2017	Yijie Li	1955404
SM 4500 F-C-97	12/20/2017			12/27/2017	Dale L. Simmons	1955404
SM 5310 B-00	12/20/2017			12/22/2017	Megan Treadwell	1955405
USGS O-2060-01	12/20/2017	12/21/2017	Hoor Shaik	12/22/2017	Julie Michelle Frank	1955410