

Benchmark EnviroAnalytical, Inc.
1001 Corporate Ave, Suite 102
North Port, FL 34289
(941) 625-3137 / (941) 240-3071 fax

Client Name: FDEP
Address: 2600 Blair Stone Rd.
Tallahassee, FL, 32399
Phone: (850) 245-8085
Email: dep-overflowmicro@dep.state.fl.us

Bill to: Frances Moore
P.O. #: Lab B16970

Project Name: RQ-2019-04-01-64

Project Run/Location: HA's - Alligator & Hudson

Laboratory Submission #:

SI9040044

Station ID	Surface Water Sample Matrix ¹ (Fresh / Marine) Sample Type	Enterococci E.coli (SM9223B Quantitray)	Laboratory Sample #
V.Coli = 250 mL ; Entero = 120 mL ; Fecal Coliform = 120 mL			
G2SD0033-Alligator Creek South	F / M / Grab	E.coli / Entero / F.coli	04/02/2019 1240 0925
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:
	F / M / Grab	E.coli / Entero / F.coli	Date & Time:

¹ Sample Type¹ is used to indicate whether the sample was a grab (G) or whether it was a composite (C).

² Sample Matrix² is used to indicate whether the sample is being discharged to drinking water (DW), groundwater (GW), surface water (SW), fresh surface water (FSW), saline surface water (SSW), soil sediment (SDMNT), or sludge (SLDG).

³ Container Type³ is used to indicate whether the container is plastic (P) or glass (G).

⁴ Sample must be refrigerated or stored in ice after collection. The maximum temperature during storage should be 6°C (41°F). Under "Preservative," list any preservatives that were added to the sample container.

INSTRUCTIONS

- Each bottle has a label identifying sample ID, preservation preservative contained in the bottle, sample type, client ID, and parameters for analysis.
- The following information should be added to each bottle label after collection with permanent black ink: date and time of collection, sampler's name or initials, and any field number or ID.
- All bottles not containing preservative may be used with appropriate sample prior to collection.
- The client is responsible for documentation of the sampling event. Please note a special sampling event on the sample custody form.

Laboratory Sample Acceptability:

pH < 2.0 BEAS Temperature: 1.9°C

1	Collector: <i>Lawrence Morgan</i>	Date: 04/02/19	Time: 1408	Received By:	Date:	Time:
2	Relinquished by: <i>W. J. B. B. B.</i>	Date: 04/02/19	Time: 1408	Received By: <i>W. J. B. B. B.</i>	Date: 4/2/19	Time: 1408
3	Relinquished by:	Date:	Time:	Received By:	Date:	Time:
4	Relinquished by:	Date:	Time:	Received By:	Date:	Time:
5	Relinquished by:	Date:	Time:	Received By:	Date:	Time: