

Group: A	Bottle Type:	P-1L	# of Bottles: 1	Preservative:	ICE	Enter Number of Bottles Sent to Lab
	ALKALINITY					
	TURBIDITY					
	W-CL-IC					
	W-COLOR					
	W-F					
	W-SO4-IC					
	W-TDS					
	W-TSS					
Group: A	Bottle Type:	BP-1L	# of Bottles: 1	Preservative:	ICE	Enter Number of Bottles Sent to Lab
	CHLSUITE-W					
Group: A	Bottle Type:	BG-500ML	# of Bottles: 1	Preservative: CH2ClCOOH and Ice		Enter Number of Bottles Sent to Lab
	W-CARB-AA		128343	2/4/23		
Group: A	Bottle Type:	BG-500ML FILL TO TOP	# of Bottles: 4	Preservative:	ICE	Enter Number of Bottles Sent to Lab
	W-CT-CI-R					
	W-PCL-TQ-R					
Group: A	Bottle Type:	P-125ML	# of Bottles: 1	Preservative: FILTER-H2SO4-ICE		Enter Number of Bottles Sent to Lab
	W-DOC			SA1308030	12/1/22	
Group: A	Bottle Type:	BG-500ML	# of Bottles: 1	Preservative:	ICE	Enter Number of Bottles Sent to Lab
	W-E8321-DI					
	W-E8321-MS					
Group: A	Bottle Type:	P-500ML	# of Bottles: 1	Preservative: HNO3		Enter Number of Bottles Sent to Lab
	W-HARD			NA1105100		
	W-ICP					
	W-ICPMS					
Group: A	Bottle Type:	P-500ML	# of Bottles: 1	Preservative: ICE And H2SO4		Enter Number of Bottles Sent to Lab
	W-NH3					
	W-NO2NO3			SA1308030		
	W-S-A-TP					
	W-TKN					
	W-TOC					
Group: A	Bottle Type:	P-125ML	# of Bottles: 1	Preservative:	FILTER-ICE	Enter Number of Bottles Sent to Lab
	W-PO4-F					

Group: A

Bottle Type:

BG-1L

of Bottles: 4

Preservative:

ICE

Enter Number of Bottles Sent to Lab

4

W-PSNP-TQ

