

Chemical Analysis Report

NE-DIST-2016-05-26-01

Florida Department of Environmental Protection
Central Laboratory
2600 Blair Stone Road
Tallahassee, FL 32399-2400
DOH Accreditation E31780

Event Description: **Reference Streams-SCI**

Request ID: **RQ-2016-05-23-12**

Customer: **NE-DIST**


Project ID: **SCIREF**

Send Reports to:
FL Dept. of Environmental Protection
8800 Baymeadows Way West
Jacksonville, FL 32256-7577
Attn: Jeremy Parrish

For additional information please contact
Timothy W. Fitzpatrick
Liang-Tsair Lin, Ph.D.
Colin Wright, Ph.D.
Thekkekalathil Chandrasekhar, PhD, QA Officer
Phone (850) 245-8085

Certified by: Liang Lin, Environmental Administrator

Date Certified: 20-JUN-2016 14:55

A handwritten signature in black ink, appearing to read "Liang Lin", is shown on a light-colored background.

Case Narrative

Unless otherwise noted, all samples included in this report were received in accordance with protocols referenced in Chapter 62-160, Florida Administrative Code (F.A.C.). Results published in this report pertain only to the samples as submitted to, and received by the laboratory. All times in this report are adjusted to the applicable Eastern Time Zone (EST or EDT).

Results for the following analytical groups are included in this report: Metals, Nutrients and Pesticides.

Scientific notation may be used in reporting very large or small values. Values reported using scientific notation will take the form of the following example: 1.3E+03, which is equivalent to 1.3×10^3 or 1300.

Unless otherwise noted, analytical values for soil and sediment samples are reported on a dry weight basis, and analytical values for waste and tissue samples are reported on a wet weight basis.

Results for TNI accredited tests met requirements established by The NELAC Institute. A double asterisk (**) is used to indicate an analyte/matrix/method for which the laboratory is not TNI accredited or where accreditation for that field of testing is not applicable.

Any significant anomalies or deviations from established protocols are documented in Non-Conformance Reports, which, where appropriate, are included within this analytical report. Additional comments related to specific analytical tests may be included as remarks following the analytical results for each sample. Such comments and remarks are for informational purposes only and are not intended to convey judgement about the usability of the reported data.

A quality control report on the performance of the test method for the submitted samples is included. Uncertainty associated with the analytical results contained in this report can be estimated from the reported quality assurance results and from published quality control acceptance limits for each analytical test. Matrix quality control results (matrix spike recoveries and matrix sample precision) pertain only to the matrix sample tested and do not necessarily reflect test method performance for other samples.

Typical matrix quality control (QC) measurements may include matrix spike recovery, matrix spike duplicate recovery, matrix spike precision and matrix sample precision. Not all matrix QC results may be available or reportable; where they are not an explanation is provided. Typical reasons for unavailable QC results include, but are not limited to, a) insufficient matrix sample to perform some or all QC measurements; b) analyte concentration in the sample replicated was too low for a meaningful measurement of precision and c) analyte concentration in the matrix sample spiked was too high (relative to the amount of analyte spiked) for a meaningful measurement of recovery. Where matrix QC results are unavailable, other method performance metrics (e.g., LCS recovery, LCS precision, surrogate recovery) may be used to assess performance of the method. Comments explaining any missing QC measurements are not intended to convey any adverse conclusions about the quality of the reported data.

Precision is reported as relative percent difference unless otherwise noted.

Quality Control codes as defined below may be used in this report to indicate results that are associated with one or more quality control elements which did not fall within established test method criteria. Such results may be qualified as estimates using a J qualifier as required by 62-160 F.A.C. Explanations are included in the report for any results that were reported as estimates for other reasons.

QC Codes used in this report may include:

- LCS – Recovery for the batch Laboratory Control Sample (LCS) was outside existing control limits;
- MS – Recovery for the batch matrix spike (MS) was outside existing control limits;
- CCV – Recovery for a continuing calibration verification (CCV) standard was outside existing control limits;
- SUR – Recovery of a surrogate (SUR) for associated analytes was outside existing control limits;
- RPD – The precision, measured as relative percent difference (RPD), of batch replicate measurements was outside existing control limits;
- RSD – The precision, measured as relative standard deviation (RSD), of batch replicate measurements was outside existing control limits;
- SMP – Sample - used precision derived from replicate analyses of a sample;

The following data qualifiers are used, where applicable, in this report as specified in 62-160 F.A.C.

- A - Value reported is the mean of two or more determinations.
- B - Results based on colony counts outside the acceptable range.
- I - The reported value is between the laboratory method detection limit and the laboratory practical quantitation limit.
- J - Estimated value and/or the analysis did not meet established quality control criteria.
- K - Actual value is known to be less than value given.
- L - Actual value is known to be greater than value given.
- N - Presumptive evidence of presence of material.
- O - Sampled, but analysis lost or not performed.
- Q - Sample held beyond normal holding time.
- T - Value reported is less than the criterion of detection.
- U - Material was analyzed for but not detected. The reported value is the method detection limit for the sample analyzed.
- V - Analyte was detected in both sample and method blank.
- X - Too few individuals to calculate SCI value.
- Y - The laboratory analysis was from an unpreserved or improperly preserved sample. The data may not be accurate.
- Z - Colonies were too numerous to count (TNTC).

Quality control information from overflow laboratories may not be included in this report. Please refer to the associated report from the overflow laboratory for additional information.

Sample Location: ATES CREEK @ 315 SHARON RD.

Collection Date/Time: 05/25/2016 10:00

Field ID: ATES CREEK @ 315 SHARON RD.

Matrix: W-SURF-FRH

Sample ID	Ref. Method	Component	Result	Code	Units	Batch ID	QC Codes
1802418	EPA 180.1 Rev. 2.0	Turbidity	1.9		NTU	P305150	
	EPA 300.0	Bromide**	0.029	I	mg Br/L	P305180	
	EPA 300.0 Rev. 2.1	Chloride	7.9		mg Cl/L	P305175	
		Sulfate	5.0		mg SO4/L	P305176	
	SM 2120 B	Color (true)	250		PCU	P305125	
	SM 2320 B-97	Total Alkalinity	0.65	U	mg CaCO3/L	P305510	
	SM 2540 C-97	TDS	64		mg/L	P305892	
	SM 2540 D-97	TSS	2	U	mg/L	P305895	
	SM 4500 F-C-97	Fluoride	0.099	I	mg F/L	P305512	
1802419	EPA 350.1 Rev. 2.0	Ammonia-N	0.017		mg N/L	P305661	
	EPA 351.2 Rev. 2.0	Kjeldahl Nitrogen	0.51		mg N/L	P305447	
	EPA 353.2 Rev. 2.0	NO2NO3-N	0.005	I	mg N/L	P305187	
	EPA 365.1 Rev. 2.0	Total-P	0.26		mg P/L	P305786	
	SM 5310 B-00	Organic Carbon	21		mg C/L	P305854	
1802420	EPA 365.1 Rev. 2.0 dissolved	O-Phosphate-P	0.22		mg P/L	P305119	
1802425	EPA 200.7 Rev. 4.4	Calcium	3.31		mg/L	P305304	
		Magnesium	0.94		mg/L	P305304	
		Potassium	0.30	U	mg/L	P305304	
		Sodium	4.2		mg/L	P305304	
1802426	EPA 8321B	Diuron	0.0040	U	ug/L	P304884	
		Imidacloprid	0.0040	U	ug/L	P304884	
		Linuron	0.0040	U	ug/L	P304884	
		Sucralose**	0.010	U	ug/L	P304884	
		Acetaminophen**	0.0040	U	ug/L	P304884	
		Carbamazepine**	4.0E-04	U	ug/L	P304884	
		Primidone**	0.0080	U	ug/L	P304884	
		Fluridone**	4.0E-04	U	ug/L	P304884	

Ref. Method and Comment:

SM 2540 D-97: The precision data is unavailable due to the low analyte concentration in the QC sample.

Quality Assurance Report

Method Blank Results

Reference Method: EPA 180.1 Rev. 2.0

Batch ID: P305150

Component	Result	Code	Units
Turbidity	0.10	U	NTU

Reference Method: EPA 200.7 Rev. 4.4

Batch ID: P305304

Component	Result	Code	Units
Calcium	0.075	U	mg/L
Magnesium	0.040	U	mg/L
Potassium	0.30	U	mg/L
Sodium	0.50	U	mg/L

Quality Assurance Report Method Blank Results

Reference Method: EPA 300.0
Batch ID: P305180

Component	Result	Code	Units
Bromide	0.016	U	mg Br/L

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305175

Component	Result	Code	Units
Chloride	0.20	U	mg Cl/L

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305176

Component	Result	Code	Units
Sulfate	0.20	U	mg SO4/L

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P305661

Component	Result	Code	Units
Ammonia-N	0.002	U	mg N/L

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P305447

Component	Result	Code	Units
Kjeldahl Nitrogen	0.080	U	mg N/L

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P305187

Component	Result	Code	Units
NO2NO3-N	0.004	U	mg N/L

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P305786

Component	Result	Code	Units
Total-P	0.002	U	mg P/L

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P305119

Component	Result	Code	Units
O-Phosphate-P	0.004	U	mg P/L

Reference Method: EPA 8321B
Batch ID: P304884

Component	Result	Code	Units
Acetaminophen	0.0040	U	ug/L
Carbamazepine	4.0E-04	U	ug/L
Diuron	0.0020	U	ug/L
Fluridone	4.0E-04	U	ug/L

Quality Assurance Report Method Blank Results

Reference Method: EPA 8321B
Batch ID: P304884

Component	Result	Code	Units
Imidacloprid	0.0040	U	ug/L
Linuron	0.0040	U	ug/L
Primidone	0.0080	U	ug/L
Sucralose	0.010	U	ug/L

Reference Method: SM 2120 B
Batch ID: P305125

Component	Result	Code	Units
Color (true)	2.5	U	PCU

Reference Method: SM 2320 B-97
Batch ID: P305510

Component	Result	Code	Units
Total Alkalinity	0.65	U	mg CaCO3/L

Reference Method: SM 2540 C-97
Batch ID: P305892

Component	Result	Code	Units
TDS	15	U	mg/L

Reference Method: SM 2540 D-97
Batch ID: P305895

Component	Result	Code	Units
TSS	2	U	mg/L

Reference Method: SM 4500 F-C-97
Batch ID: P305512

Component	Result	Code	Units
Fluoride	0.035	U	mg F/L

Reference Method: SM 5310 B-00
Batch ID: P305854

Component	Result	Code	Units
Organic Carbon	0.50	U	mg C/L

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: EPA 200.7 Rev. 4.4
Batch ID: P305304

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Calcium	106		P	85 - 115
Magnesium	105		P	85 - 115
Potassium	100		P	85 - 115

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: EPA 200.7 Rev. 4.4
Batch ID: P305304

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Sodium	104		P	85 - 115

Reference Method: EPA 300.0
Batch ID: P305180

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Bromide	101		P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305175

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Chloride	99.6		P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305176

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Sulfate	98.6		P	90 - 110

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P305661

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Ammonia-N	98.8		P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P305447

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Kjeldahl Nitrogen	98.7		P	90 - 110

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P305187

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
NO2NO3-N	100		P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P305786

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Total-P	97.8		P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P305119

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
O-Phosphate-P	104		P	90 - 110

Quality Assurance Report Laboratory Control Sample Accuracy

Reference Method: EPA 8321B
Batch ID: P304884

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Acetaminophen	120		P	40 - 140
Carbamazepine	103		P	40 - 140
Diuron	107		P	40 - 140
Fluridone	104		P	40 - 140
Imidacloprid	116		P	40 - 160
Linuron	89.2		P	40 - 140
Primidone	128		P	40 - 140
Sucralose	95.4		P	40 - 140

Reference Method: SM 2320 B-97
Batch ID: P305510

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Total Alkalinity	101		P	95 - 105

Reference Method: SM 4500 F-C-97
Batch ID: P305512

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Fluoride	98.2		P	92 - 103

Reference Method: SM 5310 B-00
Batch ID: P305854

Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
Organic Carbon	101		P	85 - 115

Quality Assurance Report Matrix Spike Accuracy

Reference Method: EPA 200.7 Rev. 4.4
Batch ID: P305304

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802409	Calcium	102	102	P/P	80 - 120
1802409	Magnesium	101	103	P/P	80 - 120
1802409	Potassium	97.1	97.9	P/P	80 - 120
1802409	Sodium	98.7	98.9	P/P	80 - 120

Reference Method: EPA 300.0
Batch ID: P305180

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1800451	Bromide	99.5	100	P/P	90 - 110
1800628	Bromide	102		P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305175

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1800451	Chloride	96.5	97.4	P/P	90 - 110

Quality Assurance Report Matrix Spike Accuracy

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305175

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1800628	Chloride	98.8		P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305176

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1800451	Sulfate	96.9	98.2	P/P	90 - 110
1800628	Sulfate	98.3		P	90 - 110

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P305661

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802367	Ammonia-N	101	102	P/P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P305447

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802340	Kjeldahl Nitrogen	102	98.6	P/P	90 - 110

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P305187

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802364	NO2NO3-N	98.5	100	P/P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P305786

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802534	Total-P	98.8	101	P/P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P305119

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802368	O-Phosphate-P	101	102	P/P	90 - 110

Reference Method: EPA 8321B
Batch ID: P304884

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1801187	Acetaminophen	72.4	76.8	P/P	40 - 140
1801187	Carbamazepine	76.0	75.2	P/P	40 - 140
1801187	Diuron	69.2	106	P/P	40 - 140
1801187	Fluridone	93.4	101	P/P	40 - 140
1801187	Imidacloprid	122	136	P/P	40 - 160
1801187	Linuron	73.6	79.6	P/P	40 - 140
1801187	Primidone	92.0	102	P/P	40 - 140
1801187	Sucralose	56.0	54.7	P/P	40 - 140

Quality Assurance Report Matrix Spike Accuracy

Reference Method: SM 4500 F-C-97
Batch ID: P305512

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802519	Fluoride	94.7	95.3	P/P	91 - 105

Reference Method: SM 5310 B-00
Batch ID: P305854

Spiked Sample	Component	% Rec.1	% Rec.2	Pass/Fail	Control Limits
1802292	Organic Carbon	101	101	P/P	85 - 115

Quality Assurance Report Precision

Reference Method: EPA 180.1 Rev. 2.0
Batch ID: P305150

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802333	Turbidity	3.53	Sample	P	0 - 20

Reference Method: EPA 200.7 Rev. 4.4
Batch ID: P305304

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802409	Calcium	0.452	Spike	P	0 - 20
1802409	Magnesium	1.82	Spike	P	0 - 20
1802409	Potassium	0.809	Spike	P	0 - 20
1802409	Sodium	0.169	Spike	P	0 - 20

Reference Method: EPA 300.0
Batch ID: P305180

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1800451	Bromide	0.658	Spike	P	0 - 20

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305175

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1800451	Chloride	0.588	Spike	P	0 - 20

Reference Method: EPA 300.0 Rev. 2.1
Batch ID: P305176

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1800451	Sulfate	1.26	Spike	P	0 - 20

Quality Assurance Report Precision

Reference Method: EPA 350.1 Rev. 2.0
Batch ID: P305661

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802367	Ammonia-N	0.894	Spike	P	0 - 20

Reference Method: EPA 351.2 Rev. 2.0
Batch ID: P305447

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802340	Kjeldahl Nitrogen	3.06	Spike	P	0 - 20

Reference Method: EPA 353.2 Rev. 2.0
Batch ID: P305187

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802364	NO2NO3-N	0.960	Spike	P	0 - 20

Reference Method: EPA 365.1 Rev. 2.0
Batch ID: P305786

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802534	Total-P	1.76	Spike	P	0 - 20

Reference Method: EPA 365.1 Rev. 2.0 dissolved
Batch ID: P305119

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802368	O-Phosphate-P	0.891	Spike	P	0 - 20

Reference Method: EPA 8321B
Batch ID: P304884

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1801187	Acetaminophen	5.90	Spike	P	0 - 30
1801187	Carbamazepine	1.06	Spike	P	0 - 30
1801187	Diuron	12.1	Spike	P	0 - 30
1801187	Fluridone	7.21	Spike	P	0 - 30
1801187	Imidacloprid	4.30	Spike	P	0 - 30
1801187	Linuron	7.83	Spike	P	0 - 30
1801187	Primidone	10.3	Spike	P	0 - 30
1801187	Sucralose	2.35	Spike	P	0 - 30

Reference Method: SM 2120 B
Batch ID: P305125

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802362	Color (true)	0.251	Sample	P	0 - 20

Quality Assurance Report Precision

Reference Method: SM 2320 B-97
Batch ID: P305510

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802519	Total Alkalinity	0.261	Sample	P	0 - 4.8

Reference Method: SM 2540 C-97
Batch ID: P305892

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802429	TDS	2.09	Sample	P	0 - 10

Reference Method: SM 4500 F-C-97
Batch ID: P305512

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802519	Fluoride	0.536	Spike	P	0 - 3.0

Reference Method: SM 5310 B-00
Batch ID: P305854

Replicated Lab Sample	Component	% RSD/RPD	Sample/Spike/LCS*	Pass/Fail	Control Limits
1802292	Organic Carbon	0.287	Spike	P	0 - 20

* Sample, spike and/or laboratory control sample precision (LCS) is reported.
Replicate spike precision may be reported when sample results are below quantifiable levels.

Quality Assurance Report Surrogates

Lab Sample ID: 1802426
Field Sample ID: ATES CREEK @ 315 SHARON RD.

Reference Method	Surrogate	% Rec.	Pass/Fail	Control Limits
EPA 8321B	Acetaminophen-d4	91.4	P	30 - 150
EPA 8321B	Acetaminophen-d4	91.4	P	30 - 150
EPA 8321B	Carbamazepine-d10	86.8	P	40 - 160
EPA 8321B	Carbamazepine-d10	86.8	P	40 - 160
EPA 8321B	Diuron-d6	85.2	P	40 - 150
EPA 8321B	Diuron-d6	85.2	P	40 - 150
EPA 8321B	Primidone-d5	108	P	40 - 150
EPA 8321B	Primidone-d5	108	P	40 - 150
EPA 8321B	Sucralose-d6	44.0	P	40 - 150
EPA 8321B	Sucralose-d6	44.0	P	40 - 150

Quality Assurance Report Calibration Verification

Quality Assurance Report Calibration Verification

Reference Method: EPA 300.0

Run ID: A69676

Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Bromide	100	101	P/P	90 - 110

Reference Method: EPA 300.0 Rev. 2.1

Run ID: A69676

Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Chloride	101	101	P/P	90 - 110
Sulfate	99.8	97.3	P/P	90 - 110

Reference Method: EPA 365.1 Rev. 2.0 dissolved

Run ID: A69721

Included Lab Sample IDs: 1802420

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
O-Phosphate-P	98.5	98.7	P/P	90 - 110

Reference Method: SM 2120 B

Run ID: A69723

Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Color (true)	99.7	100	P/P	90 - 115

Reference Method: EPA 180.1 Rev. 2.0

Run ID: A69730

Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Turbidity	100	101	P/P	95 - 105

Reference Method: EPA 8321B

Run ID: A69733

Included Lab Sample IDs: 1802426

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Acetaminophen	118	117	P/P	60 - 150
Carbamazepine	114	112	P/P	60 - 150
Diuron	116	127	P/P	60 - 150
Fluridone	108	108	P/P	60 - 150
Imidacloprid	129	138	P/P	60 - 150
Linuron	118	107	P/P	60 - 150
Primidone	129	115	P/P	60 - 150

Reference Method: EPA 353.2 Rev. 2.0

Run ID: A69735

Included Lab Sample IDs: 1802419

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
NO2NO3-N	97.0	96.5	P/P	90 - 110

Quality Assurance Report Calibration Verification

Reference Method: EPA 8321B
Run ID: A69770
Included Lab Sample IDs: 1802426

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Sucralose	80.4	84.3	P/P	70 - 150

Reference Method: EPA 200.7 Rev. 4.4
Run ID: A69794
Included Lab Sample IDs: 1802425

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Calcium	99.2	102	P/P	90 - 110
Magnesium	99.7	104	P/P	90 - 110
Potassium	97.0	99.9	P/P	90 - 110
Sodium	98.0	99.3	P/P	90 - 110

Reference Method: SM 2320 B-97
Run ID: A69825
Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Total Alkalinity	102	101	P/P	90 - 110

Reference Method: SM 4500 F-C-97
Run ID: A69825
Included Lab Sample IDs: 1802418

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Fluoride	100	101	P/P	90 - 110

Reference Method: EPA 350.1 Rev. 2.0
Run ID: A69862
Included Lab Sample IDs: 1802419

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Ammonia-N	100	101	P/P	90 - 110

Reference Method: EPA 351.2 Rev. 2.0
Run ID: A69880
Included Lab Sample IDs: 1802419

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Kjeldahl Nitrogen	97.8	95.9	P/P	90 - 110

Reference Method: SM 5310 B-00
Run ID: A69933
Included Lab Sample IDs: 1802419

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Organic Carbon	96.4	99.0	P/P	90 - 110

Quality Assurance Report Calibration Verification

Reference Method: EPA 365.1 Rev. 2.0

Run ID: A69962

Included Lab Sample IDs: 1802419

Component	% Rec.1	% Rec.2	Pass/Fail*	Control Limits
Total-P	101	102	P/P	90 - 110

* Pass/Fail determinations are made for each bracketing calibration verification check.
Control limits for initial calibration checks may be different from those for continuing checks, depending on method requirements.
Where they are different, both control limits are provided.

Quality Assurance Report Summary

Ref. Method	Analyte	LCS % Recovery	MS % Recovery	LCS	Precision SMP	MS
EPA 180.1 Rev. 2.0	Turbidity				3.53	
EPA 200.7 Rev. 4.4	Calcium	106	102	102		0.452
	Magnesium	105	101	103		1.82
	Potassium	100	97.1	97.9		0.809
	Sodium	104	98.7	98.9		0.169
EPA 300.0	Bromide	101	99.5	100	102	0.658
EPA 300.0 Rev. 2.1	Chloride	99.6	96.5	97.4	98.8	0.588
	Sulfate	98.6	96.9	98.2	98.3	1.26
EPA 350.1 Rev. 2.0	Ammonia-N	98.8	101	102		0.894
EPA 351.2 Rev. 2.0	Kjeldahl Nitrogen	98.7	102	98.6		3.06
EPA 353.2 Rev. 2.0	NO2NO3-N	100	98.5	100		0.960
EPA 365.1 Rev. 2.0	Total-P	97.8	98.8	101		1.76
EPA 365.1 Rev. 2.0 dissolved	O-Phosphate-P	104	101	102		0.891
EPA 8321B	Acetaminophen	120	72.4	76.8		5.90
	Carbamazepine	103	76.0	75.2		1.06
	Diuron	107	69.2	106		12.1
	Fluridone	104	93.4	101		7.21
	Imidacloprid	116	122	136		4.30
	Linuron	89.2	73.6	79.6		7.83
	Primidone	128	92.0	102		10.3
	Sucralose	95.4	56.0	54.7		2.35
SM 2120 B	Color (true)				0.251	
SM 2320 B-97	Total Alkalinity	101			0.261	
SM 2540 C-97	TDS				2.09	
SM 4500 F-C-97	Fluoride	98.2	94.7	95.3		0.536
SM 5310 B-00	Organic Carbon	101	101	101		0.287

Reference Method Descriptions

Method / DoH Cert #	Description	Associated Samples
EPA 180.1 Rev. 2.0 / E31780	Turbidity in aqueous matrices	1802418
EPA 200.7 Rev. 4.4 / E31780	Total Recoverable Metals analysis using ICP emission spectroscopy for aqueous samples supporting Clean Water Act Projects	1802425
EPA 300.0 / E31780	Bromide in aqueous matrices	1802418
EPA 300.0 Rev. 2.1 / E31780	Chloride in aqueous matrices	1802418
EPA 300.0 Rev. 2.1 / E31780	Sulfate in aqueous matrices	1802418

Reference Method Descriptions

Method / DoH Cert #	Description	<u>Associated Samples</u>
EPA 350.1 Rev. 2.0 / E31780	Ammonia in aqueous matrices as mg N/L	1802419
EPA 351.2 Rev. 2.0 / E31780	Total Kjeldahl Nitrogen in aqueous matrices	1802419
EPA 353.2 Rev. 2.0 / E31780	Nitrite/Nitrate in aqueous matrices as mg N/L	1802419
EPA 365.1 Rev. 2.0 / E31780	Total Phosphorus in aqueous matrices as mg P/L	1802419
EPA 365.1 Rev. 2.0 dissolved / E31780	Ortho-phosphate, dissolved, in filtered, aqueous matrices as mg P/L	1802420
EPA 8321B / E31780	Analysis of sucralose in water matrices by HPLC/MS/MS	1802426
EPA 8321B / E31780	Personal Care Products and pesticides in water matrices by HPLC/MS/MS	1802426
SM 2120 B / E31780	True Color measured at 450 nm	1802418
SM 2320 B-97 / E31780	Alkalinity in aqueous matrices as mg CaCO ₃ /L	1802418
SM 2540 C-97 / E31780	Total Dissolved Solids in aqueous matrices	1802418
SM 2540 D-97 / E31780	Total Suspended Solids in aqueous matrices	1802418
SM 4500 F-C-97 / E31780	Fluoride in aqueous matrices, without distillation (not valid for NPDES monitoring)	1802418
SM 5310 B-00 / E31780	Nonpurgeable Organic Carbon in aqueous matrices	1802419

Preparation and Analysis Log

Ref. Method	Received Date	Prep Date/Time	Prepared By	Analysis Date/Time	Analyzed By	Associated Samples
EPA 180.1 Rev. 2.0	05/26/2016			05/26/2016 16:42	Blanca Fach	1802418
EPA 200.7 Rev. 4.4	05/26/2016	05/31/2016	Elliott D. Healy	06/01/2016	Joshua Ayres	1802425
EPA 300.0	05/26/2016			05/26/2016	Ping Hua	1802418
EPA 300.0 Rev. 2.1	05/26/2016			05/26/2016	Ping Hua	1802418
EPA 350.1 Rev. 2.0	05/26/2016			06/06/2016	Diana M. Turner	1802419
EPA 351.2 Rev. 2.0	05/26/2016	06/02/2016	Sofia Elnemr	06/06/2016	Christopher Armour	1802419
EPA 353.2 Rev. 2.0	05/26/2016			05/27/2016	Peter D. Kaus	1802419
EPA 365.1 Rev. 2.0	05/26/2016	06/08/2016	Vijayalakshmi Reddy	06/09/2016	Peter D. Kaus	1802419
EPA 365.1 Rev. 2.0 dissolved	05/26/2016			05/26/2016 15:25	Julia Storbeck	1802420
EPA 8321B	05/26/2016	05/26/2016	Hoor Shaik	05/27/2016	Juyoung Kim	1802426
	05/26/2016	05/26/2016	Hoor Shaik	05/31/2016	Juyoung Kim	1802426
SM 2120 B	05/26/2016			05/26/2016 17:18	Chrisoulla Rakowski	1802418
SM 2320 B-97	05/26/2016			06/03/2016	Dale L. Simmons	1802418
SM 2540 C-97	05/26/2016			05/31/2016	Blanca Fach	1802418
SM 2540 D-97	05/26/2016			05/31/2016	Blanca Fach	1802418
SM 4500 F-C-97	05/26/2016			06/03/2016	Dale L. Simmons	1802418
SM 5310 B-00	05/26/2016			06/08/2016	Diana M. Turner	1802419