

**FLORIDA DEPARTMENT OF ENVIRONMENTAL PROTECTION**

STRATEGIC MONITORING PLAN Field Sheet

July 2022

Meter
Passed End
of Day CCV
Yes / No

WIN Station Name:		TWIN LAKE				Sample Date:		4/10/2023	
WIN / Field ID:	G3NW0215	ORG ID:	21FLPNS	ENR:		Latitude:	30.395162		
County:	OKALOOSA	WBID:	906A	LAKE	3F	Longitude:	-86.47528		
Receiving Body of Water:		Coleman Lake				RQ-	2023-04-10-36		

Sampling Team Members	WQ Measurements	Sample Collection	Documentation	Preservation	Preservation Check
Frank Butera	x		x		
Hannah Henning		x		x	x

Sampling Equipment			
x	Sample Container		Carboy
	Dipper		Van Dorn
	Syringe		Pole
Syringe Lot#		Filter Lot#	
Secchi Equipment ID:		BFA	
Collection Method			
	Wading	x	Canoe/Kayak
	From Shore		Electric Motor
			Gasoline Motor
Depth Measured With		Meter	

Water Level:	Normal
Flow:	NA
Matrix:	Fresh
Tide Stage:	
High Tide:	
Low Tide:	
Photos:	NA
Meter ID:	154179 (ProDSS)
Depth Gauge ID:	depth_gauge

Time Zone	Time	Total Depth (m)	Sample Depth (m)	DO (mg/L)	DO (%SAT)	Temp. (°C)	pH	Secchi Disk (m)	Sp COND (µmhos/cm)	Salinity (PPTth)
CEN	1015	0.47	0.24	2.03	21.9	19.2	6.26		257.6	0.12

Biological Samples Collected: None

SBIO Visit ID: HA ID: RPS ID: Macrophyte ID: LIMS Sample ID Entered Into The Macrophyte Module:

Parameter Suite	Parameter Methods (LIMS Test Codes)	Preservation				
		Type	Lot#	Expiration Date	pH Check	
Microbiology						
Chlorophyll	CHLSUITE-W	ICE				
Physical/Aggregate Nutrients	ALKALINITY, TURBIDITY, W-F, W-TSS W-CL-IC, W-COLOR, W-SO4-IC, W-TDS	ICE				
Filtered Nutrients			Syringe Lot#; Filter Lot#:			
Acid Preserved Nutrients	W-NH3 / W-NO2NO3 / W-TKN / W-TOC / W-S-A-TP	ICE	H2SO4*	SA2229040	11/04/2023	H2SO4 pH<2
Metals						
Organics (LC)						
Macroinvertebrates						
Organics (GC)			*2.0 ml of H2SO4 Used Unless Noted Otherwise			
Molecular						
Periphyton						
Phytoplankton						
BOD						
DOC						
Other(s)			SyringeLot#		FilterLot#	
Field Comments:	Veg. too thick to get secchi measurement					
Comments on COC:						
Relinquish Date:	04/10/2023	Signature:				

LIMS EVENT ID: WIN Import ID:

Enter Field Parameters, Verify Station ID In LIMS DMT: QC Station ID, Field Parameters In LIMS DMT:

Save Field Sheets on Network: Migrate Data to WIN: