Temperature Log for Solids

DEP SOPs NU-047 (TDS), NU-051 (TSS) and NU-099 (TS)

Oven ID #: _____

Filename: _____ Thermometer Calibration Date: _____

Thermometer S/N: _____ Thermometer Expiration Date: _____

For Total Dissolved Solids (TDS) analysis, the oven must be 180 ± 2 °C.
For Total Solids (TS) and Total Suspended Solids (TSS) analysis, the oven must be 103 - 105 °C.

Date	Analyst Initials	START Time	STOP Time	MIN temp (°C)	MAX temp (°C)	Oven in compliance? (Y/N)	Corrective action? (Y/N)	Comments