Daily Balance Use Record of Mettler Toledo XA105DU
(Weights Bracketting Measured Weights)

Check Weight IDs
BIO-R-WEIGHTS-03 = R03
BIO-R-WEIGHTS-04 = R04
BIO-R-WEIGHTS-01 = R05
BIO-S-WEIGHTS-01 = S01

STATE of FLORIDA
DEPARTMENT of ENVIRONMENTAL PROTECTION
BIOLOGY PROGRAM

Mettler Toledo XA105DU, Serial Number B314195659

Date	Material Weight or Weight Range	Lower Weight Value	Lower Weight Measured Value	Upper Weight Value	Upper Weight Measured Value	Check Weight ID	Weighed by: (Initials)	Comments

NOTE: Record g for grams; mg for milligrams
1. One lower value and one upper value check-weight can be used to bracket an event consisting of numerous measurements.
Additional check-weight measurements are required (per day) if more than two hours elapse between measurements.

Weight Tolerance - ≤ 5g = ± 0.2mg **:** > 5g ≤ 100g = ± 0.3mg