Pseudokirchneriella subcapitata for *Ceriodaphnia dubia*

Mechanical pipette used for counts _____

Batch #	Date Prepared	Count	Count Date	Initials

Cell count is raw count times the dilution factor for the Coulter counter.
Dilute or concentrate batch as required to obain cell count between 3.9x10E7 and 4.6x10E7 per mL.

Pseudokirchneriella subcapitata for *Ceriodaphnia dubia*

Cell count is raw count times the dilution factor for the Coulter counter.
Dilute or concentrate batch as required to obain cell count between 3.9x10E7 and 4.6x10E7 per mL.