FDEP Biology Laboratory-Toxicity Section
Chronic SRT Bench Sheet

Start Date: _____ Time: _____
End Date: _____ Time: _____

Organism Age:
_____
Organism Batch #:
_____

Test Type
☐ Screen
☐ Definitive
☐ Static
☐ Flow-through

Test Species: *Ceriodaphnia dubia*
(SOP: TA-07.05, Method:2002.0)

Toxicant	Batch/LIMS ID #	Days Used

Diluent	Batch #/LIMS ID #	Days Used

Test Randomizer #:	
Replicate	Culture vial
1	
2	
3	
4	
5	
6	
7	
8	
9	
10	

Food	Days Used
YCT Batch # _____	
YCT Batch # _____	
P. subcapitata Batch # _____	
P. subcapitata Batch # _____	

C. dubia fed 0.2 ml of *YCT and and *P. subcapitata* in a 1:1 ratio. YCT= yeast, cereal leaves and trout chow

Chamber Size: 30 mL
Test Volume: 20 mL

Light Intensity: 50-100 ft candles (recommended)
Photoperiod:16 hr light, 8 hr dark (recommended)

Concentration (g/L)	Mean neonates per female	% survival
Control		

Test Acceptability Criteria:		
Mean control reproduction ≥ 15 neonates per female:	Yes	No
3rd brood in 60% of surviving females:	Yes	No
Control Survival ≥ 80%:	Yes	No
Temperature Range ≤ 3° C:	Yes	No

Stats by/Date:_____
Stats verified by/Date:_____

Results	
Control % CV:	
Dunnetts PMSD (%):	
ICP25 (g/L):	
Reproduction NOEC (g/L):	
Reproduction LOEC (g/L):	
Survival NOEC (g/L):	
Survival LOEC (g/L):	

Reviewer:_____

Comments:

FDEP Biology Laboratory-Toxicity
SRT Water Quality Information Sheet

Method: EPA-821-R-02-012

Water Quality Parameters	MHWCD LIMS ID # _____ Batch # _____	MHWCD LIMS ID # _____ Batch # _____	Method	SOP	Measured By:	Transcription Checked by:
Lab Total Residual Cl_2 (mg/L)			HACH	TA-06.22		
Alkalinity (mg/L as $CaCO_3$)			HACH	TA-06.08		
Hardness (mg/L as $CaCO_3$)			HACH	TA-06.08		
Total Ammonia (mg/L as N)			HACH	TA-06.35 or TA-06.33		
Conductivity (µmhos)			FISCHER	TA-06.36		
Days Used						

Temperature Range (°C)	
Incubator #_____	
Min._____	**Max._____**
Test	
Min._____	**Max._____**
Room B246	
Min._____	**Max._____**

Comments:

FDEP Biology Laboratory - Toxicity Section

Chronic SRT Bench Sheet

Ceriodaphnia dubia Survival and Reproduction

Conc. (g/L)	Day	Replicate Number 1	2	3	4	5	6	7	8	9	10	Total neonates	Surviving original organisms	Neonates per female	Counts by:	Conc. prep by:	Fed by/Time fed:
Control	0*																@
	1																@
	2																@
	3																@
	4																@
	5																@
	6																@
	7																@
	8																@
	Total																@

Conc. (g/L)	Day	Replicate Number 1	2	3	4	5	6	7	8	9	10	Total neonates	Surviving original organisms	Neonates per female
	0*													
	1													
	2													
	3													
	4													
	5													
	6													
	7													
	8													
	Total													

Conc. = Concentration

s=split brood

(#) = 4th brood; not counted

*Day 0 is test initiation day

X = Dead (original organism)

Pipet used for feeding = P1

Loaded By: _____

Comments:

FDEP Biology Laboratory - Toxicity Section

Chronic SRT Bench Sheet

Conc. = Concentration

Ceriodaphnia dubia Survival and Reproduction

Conc. (g/L)	Day	Replicate Number 1	2	3	4	5	6	7	8	9	10	Total neonates	Surviving original organisms	Neonates per female
	0*													
	1													
	2													
	3													
	4													
	5													
	6													
	7													
	8													
	Total													

Conc. (g/L)	Day	Replicate Number 1	2	3	4	5	6	7	8	9	10	Total neonates	Surviving original organisms	Neonates per female
	0*													
	1													
	2													
	3													
	4													
	5													
	6													
	7													
	8													
	Total													

Comments:

s= split brood
(#) =4th brood; not counted
*Day 0 is test initiation day

X = Dead (original organism)

FDEP Biology Laboratory - Toxicity Section

Chronic SRT Bench Sheet

Ceriodaphnia dubia Survival and Reproduction

Conc. = Concentration

Conc. (g/L)	Day	Replicate Number										Total neonates	Surviving original organisms	Neonates per female
		1	2	3	4	5	6	7	8	9	10			
	0*													
	1													
	2													
	3													
	4													
	5													
	6													
	7													
	8													
	Total													

Conc. (g/L)	Day	Replicate Number										Total neonates	Surviving original organisms	Neonates per female
		1	2	3	4	5	6	7	8	9	10			
	0*													
	1													
	2													
	3													
	4													
	5													
	6													
	7													
	8													
	Total													

Comments:

s=split brood
(#) =4th brood; not counted
*Day 0 is test initiation day
X = Dead (original organism)

FDEP Biology Laboratory - Toxicity Section

Chronic SRT Bench Sheet

Ceriodaphnia dubia Parameters

Replicate M = "Metered" vial, a surrogate vial for parameters

Test Day ↓	Conc. g/L	New Sample Replicate	New Sample pH (su)	New Sample Temp (°C)	New Sample DO (mg/L)	New Sample Conductivity (µmhos/cm)	Old Sample Replicate	Old Sample pH (su)	Old Sample Temp (°C)	Old Sample DO (mg/L)	Old Sample Conductivity (µmhos/cm)
Day 0	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 1	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 2	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 3	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										

Conc. = Concentration, Temp = Temperature, DO = Dissolved oxygen

Comments:

FDEP Biology Laboratory - Toxicity Section
Chronic SRT Bench Sheet

Ceriodaphnia dubia Parameters

Replicate M = "Metered" vial, a surrogate vial for parameters

Test Day ↓	Conc. g/L	New Sample Replicate	New Sample pH (su)	New Sample Temp (°C)	New Sample DO (mg/L)	New Sample Conductivity (µmhos/cm)	Old Sample Replicate	Old Sample pH (su)	Old Sample Temp (°C)	Old Sample DO (mg/L)	Old Sample Conductivity (µmhos/cm)
Day 4	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 5	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 6	Control	M									
		M									
		M									
		M									
		M									
		M									
	Measured by:										
	Recorded by:										
Day 7	Control										
	Measured by:										
	Recorded by:										

Conc. = Concentration, Temp = Temperature, DO = Dissolved oxygen

Comments:

FDEP Biology Laboratory
Chronic SRT Bench Sheet
Ceriodaphnia dubia Parameters

Test ID: Dates:

METER / DEVICE	Instrument Code
DO Meter	
pH Meter	
Conductivity Meter	
Ammonia Kit	
Traceable Thermometer	
Chlorine Kit	
Alkalinity Kit	
Hardness Kit	

Entered by and date: